

Parata PASS Ware™ Packaging Software



The software that drives Parata PASS packagers exceeds expectations for flexibility and simplicity.

Parata PASS packaging makes it easy for your patients to take the right medication at the right time, every time. Parata PASS Ware makes it easy to prepare those medications safely, efficiently and accurately.

Enjoy a flexible, easy-to-use solution.

PASS Ware was developed in close partnership with pharmacists to include the features and functionality you need. PASS Ware's intuitive, attractive interface makes it easy to:

- Support cycles from 24 hours to 30 days
- Customize pouch design by patient or facility
- Add header or footer bags to separate orders
- Print reminder bags for non-oral solids or additional information
- Include optional barcode for bedside scanning
- Use a wide range of standard reports

Improve safety from script to patient.

Parata PASS Ware supports safety and accuracy from dispensing to administration. To add a new drug to Parata PASS, scan a stock bottle and PASS Ware accesses a library containing names, strengths and drug information for more than 14,000 unique NDCs. Our smart canisters communicate with PASS Ware so the PASS packager recognizes and dispenses from the correct canister, regardless of its location in the unit. During replenishment, the barcode scanner confirms a match between the drug and canister.

In long-term care communities, barcodes on PASS pouches help ensure the right patient receives the right medication at the right time.

Expect easy integration.

Parata interfaces with more than 100 pharmacy host systems, making Parata PASS an easy addition to your pharmacy. And, it's part of a complete adherence packaging solution that also includes Parata PASS Safe Loader™ and Parata Check PASS™.

Easily customize orders by facility or patient.

PASS Ware lets you filter by facility, wing, patient, medication, geographic location or packaging date to create a batch. Then, use the drag-and-drop design tool to customize each pouch with user-friendly hours of administration images or barcodes for bedside scanning.

Enjoy PASS Ware's robust and flexible reporting.

We worked closely with pharmacists to build reports that meet real-world requirements. To guide replenishment or to view drug usage and inventory, run a standard PASS Ware report. These can also be sorted, printed and exported. Our Technical Assistance Center is available to help you create reusable custom reports — just give us a call.

Take growth in stride.

Parata PASS Ware extends the capacity of PASS units beyond the standard configuration, so you can add canisters as your business or formulary grows. Simply store additional canisters outside the unit and swap them in as needed. This also reduces the need for special tablet system trays, increasing your accuracy and efficiency.

Keep pace with the changing industry.

PASS Ware is owned and developed by Parata. Our software engineers continuously innovate in response to customer feedback, new pharmacy business models and the rapidly changing health care industry.

Parata PASS Ware™ Packaging Software

THE PASS WARE SOFTWARE PACKAGE ALSO INCLUDES THE FOLLOWING HARDWARE:

- Standard PC tower
- 22-inch widescreen monitor
- Keyboard, mouse and barcode scanner

ELECTRICAL REQUIREMENTS

- Connection to an uninterruptable power supply (UPS)
- 120V, 15A, 60Hz electrical outlets



CTPASSW-SPEC-0714-1k

©2014 Parata Systems, LLC. All rights reserved.

Parata, Parata PASS Ware, Parata PASS Safe Loader and Parata Check PASS are trademarks or registered trademarks of Parata Systems, LLC.