



Parata Perl® Pouch Inspector with Cut & Roll

Key Features

Fast: Parata Perl Cut & Roll inspects and neatly bundles 1.6 pouches per second.

Efficient: With two powerful, high-speed cameras, Perl images and identifies the contents of each pouch.

Smart: Perl's intelligent, web-based inspection application—accessible from any computer or mobile device in your network—identifies existing drugs and learns new ones in only seconds.

Standard Configurations

Parata Perl Pouch Inspector with Cut & Roll

Cut & Roll Capacity

350 pouches for unit dose 55mm dispensing

275 pouches for a typical 75mm multidose pharmacy

200 pouches for fuller multi dose 95mm dispensing

Cut & Roll Inspection Speed

Up to 2.5 pouches per second

Cut & Roll Function

Cuts and tapes by patient, facility, or hour of administration

Included with Your Purchase

- Parata Perl Pouch Inspector device
 - Parata Perl Pouch Inspector web application license (2)
 - Keyboard with integrated track pad (1)
 - Correction camera (1)
 - Barcode scanner (1)
 - White/clear paper, Part #H104 (2 Rolls)
 - Bins (2) 1 Blue for alert rolls, 1 green for OK rolls
 - 4 Spools and 1 winder based on configuration choice
- Optional Complements
- Additional floor or desktop winders
 - Additional Parata Perl Pouch Inspector web application licenses

Network Requirements


Ethernet RJ-45 (Cat VI)

Power Requirements

120VAC, 15A, 60Hz, single phase, standard outlet

Parata Perl® Pouch Inspector with Cut & Roll

Perl with Cut & Roll on workflow table.



Cut, roll, and secure strips by patient, hour of administration, or facility.


You can teach the software new pills in just seconds.



Compact design needs less than 3½ square feet of counterspace.

Up to 5 concurrent users can access the web-based inspection application from any computer on the network.

Cut & Roll can bundle up to 150 medication pouches for a single patient.



Flagged pouches are routed to a separate bin for fast, easy retrieval.

Two high-speed cameras capture images so clear only a small percentage of pouches need a closer look.

No waiting. No buttons. Just swipe a strip to get started.

Each spool holds up to 1,200-2,000 medication pouches.

Dimensions

	Width (in)	Height (in)	Depth (in)	Weight (lb)
Parata Perl with Cut & Roll and 19" Spool	66.2	62.5	27.5	249
Parata Perl with Cut & Roll and 27" Spool	70.2	62.5	27.5	250